

# DWA Active Plus CHG

**FAST ACTING, ALCOHOLIC INSTRUMENT, HAND AND SURFACE DISINFECTANT PRODUCTS**

- ✓ **FAST ACTING**
- ✓ **CLEANS AND DISINFECTS IN ONE PHASE**
- ✓ **UNIVERSALLY APPLICABLE, SUITABLE FOR INSTRUMENTS, HANDS AND SURFACES**
- ✓ **WITH 2% CHLORHEXIDINE - DIGLUCONATE ACTIVE INGREDIENT**
- ✓ **CONTAINS 70% ALCOHOL**
- ✓ **ALDEHYDE-FREE**

**CE 2409**



## PRODUCT DESCRIPTION

A DWA Active Plus CHG is a wide antimicrobial spectrum instrument disinfectant bactericide, fungicide, virucide wipe cloth, with excellent cleaning and disinfecting effect. It is a product intended for the one-phase disinfection of noninvasive medical devices. The product is fragrance free, it is generally applicable on all alcohol resistant medical device surfaces that do not come into contact with mucus membranes.

## RECOMMENDED AREAS OF APPLICATION

- Disinfection of simple glass instruments and imaging diagnostic devices applied in medical care, e.g. thermometers, ultrasound heads, X-ray equipment that come into contact with skin
- Disinfection of incubators
- Disinfection of medical equipment surfaces, including glass and metal
- Applicable for hygienic hand disinfection, as well as for disinfecting alcohol resistant instruments and work surfaces in healthcare, public institutions, dentist offices, laboratories.

## SPECIAL FEATURES

The concentration of the disinfectant liquid's biocide ingredients are adjusted so it not only achieves the prescribed microbe-count reduction for the targeted contact times, but it is also not harmful to the surface type chemically or physically. The microbial effect range and contact times are certified by numerous accredited laboratory test results. The wet disinfectant wiping cloth contains the amount of impregnating disinfectant substance, which spreads evenly on the surface for disinfection, when it is applied to the surface and properly pressed, thus the biocide liquid directly contacts the microbial contaminants. The effectiveness of hand disinfection is increased by the added chlorhexidine - digluconate active ingredient.

## EFFECT RANGE

TYPE	STANDARD	CONTACT TIME
Bacteria (bactericid)		
Enterococcus hirae	EN 1276:2010	30 sec
Escherichia coli (E. coli)		
Pseudomonas aeruginosa		
Staphylococcus aureus		
Enterococcus hirae	EN 13727:2012	
Pseudomonas aeruginosa		
Staphylococcus aureus		
Pseudomonas aeruginosa	EN 14561:2006	
Enterococcus hirae		
Staphylococcus aureus		
Staphylococcus aureus	EN 16615:2015	
Enterococcus hirae		
Pseudomonas aeruginosa		
Fungus (fungicid, yeasticid)		
Candida albicans	EN 1650:2008	30 sec
Aspergillus brasiliensis		
Candida albicans	EN 13624:2014	
Candida albicans	EN 14562:2006	
Candida albicans	EN 16615:2015	
Virus (virucid)		
Adenovirus Type 5	EN 14476:2013	1 min
Murine Norovirus		
Mycrobacteria (tuberkulocid)		
Mycobacterium Terrae	EN 14348:2005	1 min

## LOGISTICS

NAME	PACKAGING	PRODUCT/BATCH	CODE
DWA Active Plus CHG liquid	0,5 l		
	1 l		
	5 l		
	10 l		
	20 l		
DWA Active Plus CHG bucket	600 sheet		
DWA Active Plus CHG cloth box	100 sheet		
DWA DWA Active Plus CHG cloth box refill	100 sheet		

## MATERIAL COMPATIBILITY

	SURFACE MATERIAL	COMPATIBILITY
METAL	Aluminum alloy 6061	Compatible
	Copper alloy 260	
	Chromium plated metal	
	Brass alloy 110	
	Nickel	
	Nickel plated metal	
	Stainless steel 304	
	Stainless steel 316	
	Steel SAE 1008/1010	
	Titanium	
	Zinc	
POLYMERS	Acryl	Compatible
	Acrylonitrile Butadiene Styron (ABS) Cycloid M MG47	
	Polybutylene terephthalate (PBT) Valox M 215 HPR	
	Polycarbonate	
	Polyester	
	Polyethylene	
	Polypropylene	
	Polystyrene	
RUBBER	PVC	Compatible
	Santoprene™	
CERAMICS	Silicone rubber	Compatible
	Glazed ceramic tile	
GLASS	Conventional glass	Compatible

## MATERIAL COMPATIBILITY

MATERIAL	COMPATIBLE
Cast iron/soft iron	Excellent
304 stainless steel	
Buna	
Fluorocarbon	
Fluoroelastomer (FKM)	
Nitrile (T5)	
Polychloroprene	
PTFE	

## INGREDIENTS

Biocide active ingredients: 20 g/kg Chlorhexidine - digluconate, 700 g/kg propane-2

## TYPE

Medical category II. a. instrument disinfectant, PT1, PT2 category biocide product.

## HAZARDS

Highly flammable liquid and steam. Severe eye irritant. May cause drowsiness or dizziness. Harmful to aquatic environment, causes permanent damage. It is best to wipe severely contaminated surfaces before disinfection. Do not use on hot surfaces!

Do not mix with other cleaning and disinfectant agents! Only use on electronic devices after complete drying! Avoid contact with open wounds.

**Always take care of safety when applying biocides! Always read the user manual and the label before use!**

